DARK 5

TECHNICAL SHEET

POWDER FREE NITRILE GLOVES
MEDICAL GRADE

 # INS-GL-DRK-100

The DARK 5 powder-free nitrile gloves offer superior quality protection for medical applications.

The 5.0 mil black nitrile gloves are a wise choice for any professional in search of reliable and durable protection.

Made from high-quality materials, these gloves offer superior resistance to tearing, punctures, and chemicals, making them ideal for use in high-risk environments.

These nitrile medical gloves are suitable for extended use. They provide excellent protection against bacteria and better resistance to a wide variety of chemicals, including petroleum-based products, oils, greases, acids, hydrocarbons, and more. Whether you work in the food, medical, automotive or any other industry requiring adequate protection, these powder-free gloves are designed for a comfortable fit and maximum dexterity while providing reliable protection.



FEATURES

- 100 gloves/box
- For single use only
- Rolled cuff keeps the glove snug at the wrist and prevents liquid from entering.
- Latex-free
- Hypoallergenic
- Non-sterile
- Ambidextrous
- Each batch of gloves is air-tested to detect pinholes

- Compliance with ASTM D6319, EN 455 Parts 1, 2, 3, 4
- ISO 9001 and ISO 13485 quality assurance
- Health Canada Approved

INDUSTRIES WHERE MOST USED

MEDICAL
FOOD SERVICE
JANITORIAL
AUTOMOTIVE
SAFETY
BEAUTY & COSMECTIC
PHARMACEUTICAL
LABORATORY

DIMENSIONS		
Code + Size	**Palm width (mm)**	**Length (mm)**
INS-GL-DRK-100S	80-89	Min.240
INS-GL-DRK-100M	90-99	Min.240
INS-GL-DRK-100L	100-109	Min.240
INS-GL-DRK-100XL	110-119	Min.240

THICKNESS	
Measured part	**Single wall (mm)**
Finger	min. 0.13
Palm	min. 0.11
Cuff	min. 0.08

PHYSICAL PROPERTIES		
Property	**Unaged**	**Aged**
TENSILE STRENGTH	min.14 Mpa	min.14 Mpa
ELONGATION	min.500 %	min.400 %
FORCE AT BREAK	min.6 N	min.6 N

INSTRUCTION FOR USE

Before putting on gloves hands must be clean and dry. Gloves do not replace the need for hand hygiene. Hands must be washed every time gloves are removed to avoid contamination. Gloves are for single use. Reuse of gloves has been associated with transmission of antibiotic-resistant organisms. Change gloves if moving from a contaminated area to a non-contaminated area. Be sure to wwear gloves that fit properly, gloves should fit snugly around wrists and hands. Different sizes are available.

100 pieces per box, 10 boxes per shipping carton. Store in a cool, dry place. Avoid excessive heat (40°C, 104°F). Shield open boxes from direct sun and fluorescent lighting. To order, please refer to a specific size and mention the product code.

INSTINCT®